




erbe

DE	Argonplasma-Applikatoren	3
EN	Rigid Argon Plasma Coagulation Applicators	17
USA	Rigid Argon Plasma Coagulation Applicators	31
FR	Applicateurs de coagulation par plasma d'argon	45
ES	Aplicadores de coagulación con plasma de argón	59
IT	Applicatori a plasma d'argon	73
PT	Aplicadores de plasma de argônio	87
EL	Εφαρμογείς πλάσματος αργού	101
NL	Argonplasma-applicatoren	117
DA	Argonplasma-applikatorer	131
SV	Argonplasma-applikatorer	145
FI	Argonplasma-aplikaattorit	159
PL	Aplikatory plazmy argonowej	173
CS	Aplikátory argonové plazmy	187
HU	Argonplazma-aplikátorok	201
RU	Аппликаторы для аргонплазменной коагуляции	215
TR	Argon plazma aplikatörleri	231
ZH	氩气电极	245
KO	아르곤 플라스마 어플리케이터	257

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DE

Argonplasma-Applikatoren

20132-250, 20132-251, 20132-252, 20132-253, 20132-254,
20132-255, 20132-256

STERILEEO

Inhalt

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NOTES ON USE

EN

Rigid Argon Plasma Coagulation Applicators

20132-250, 20132-251, 20132-252, 20132-253, 20132-254,
20132-255, 20132-256

STERILE EO




Table of Contents

IMPORTANT!	19
1 Product description	19
2 Intended use	20
3 Normal use	20
4 Maximum electrical capacity	20
5 Safety instructions	21
6 How to use	25
6.1 Checking the expiration date	25
6.2 Checking and opening the packaging	25
6.3 Checking the product	25
6.4 Connecting the product and selecting unit settings ..	25
6.5 In the case of endoscopic application: Selecting the trocar	28
6.6 Using the product	28
7 Required accessories	29
8 Disposal	29
9 Symbols	30


IMPORTANT!

Please read all information carefully.

These notes on use do not replace the user manual of the APC unit used! Read the user manual of the APC unit and ask Erbe or your distributor in case of doubt!

1 Product description

- 1 Yellow button: Activates CUT
- 2 Blue button: Activates COAG
- 3 Light gray ReMode button: Toggles between two setting options of a program (only possible with the following unit combinations: APC 2 with VIO 300 D/VIO 200 D, APC 3 with VIO 3)
- 4 Rotary wheel: for freely adjustable movement of the HF electrode
- 5 Shaft tube: for turning spatula electrodes
- 6 APC connector: for connection to the APC unit


2 Intended use

20132-256: The Rigid Argon Plasma Coagulation Applicator is intended for monopolar coagulation and devitalization of tissue using argon plasma. It is used in open surgical procedures.

20132-250, 20132-251, 20132-252, 20132-253, 20132-254 und 20132-255: The Rigid Argon Plasma Coagulation Applicators are intended for monopolar cutting and coagulation of tissue with and without argon gas or for coagulation and devitalization using argon plasma:

- used in open surgery 20132-250, 20132-251, 20132-252, 20132-253.
- used in laparoscopic surgery 20132-254, 20132-255.


3 Normal use

The argon plasma coagulation applicator is intended for use with an APC 2 / APC 3 in combination with an Erbe VIO electrosurgical unit.

The APC 2 must be equipped with an APC socket with finger switch activation. The APC 3 is already equipped with a suitable socket ex works.


The product is intended for single use.

4 Maximum electrical capacity

The maximum electrical capacity of these products is:

- 4300 V_p

5 Safety instructions

WARNING!

5.1 Sterile product for single use

Erbe recommends that the product should not be processed after use. Processing can affect the material properties and/or function of the product such that use in accordance with the intended use is no longer possible. Also, it cannot be ensured that the product will be sterile after processing. Risks include:

- Components/connections that are damaged by high temperatures (e.g. in the autoclave).
- Components (e.g., filter, extremely small lumen, etc.) that may prevent internal cleaning/disinfection.

5.2 User requirements

This product may be used only by trained medical staff, who have been shown how to use it according to the notes on use.

Only use the product if you are familiar with its handling in conjunction with the respective APC and electrosurgical unit!

5.3 General information about protection of the user, patient, and product

Do not use in the presence of combustible or explosive materials!

Never lay this product on the patient or in his/her direct vicinity.

In the case of a long activation time and/or activations which follow in quick succession (with no cooling time) the distal instrument tip can become extremely hot. This is not a material fault, rather a consequence of heavy-duty use.

Protect the product from all physical damage! Do not throw! Do not apply force of any kind!

Do not kink the tube or wrap it round the product to which it belongs!

If damaged, do not use this product!

Checks other than those described in these notes on use do not need to be carried out.



5.4 Danger of fire in bodily cavities

Before using the APC always ensure that there are no combustible or fire-supporting gases in the lumen of the organ involved.

5.5 Gas embolisms, gas emphysema


To prevent gas embolisms, the argon flow rate must not be set high enough that argon is blown into open vessels. To prevent gas embolisms or emphysema, never point the distal end of APC applicators/APC probes directly toward open vessels or press against tissue.

5.6 Danger of unintentional damage to nerve structures

Bear in mind that during an APC application in the direct vicinity of nerve structures stimulation may be induced by electric currents.

In the case of nerves which require protection it is advisable to keep the APC activation time short and adapt maximum power output to the required indication.

5.7 Risks in the case of endoscopic application

Risk of capacitive coupling: HF current can be transferred from the activated instrument to the trocar, even if the two are not touching (known as "capacitive coupling"). Risk of burns to the patient and/or operating theater staff. Do not use hybrid trocars (e.g. a combination of metal and plastic). Only use trocars designed for use with electrosurgical devices.

Unwanted stimulation of nerves and muscles due to low-frequency currents. Dangerous summation of patient leakage currents from several power sources: It is not possible to entirely avoid low-frequency currents in electrosurgery, but they can be limited by applying technical measures to equipment and accessories. Limit types show how wide-ranging these measures are, e.g., type B, BF or CF as per EN 60601. Type CF offers the highest level of protection against low-frequency currents. However, please note the following: If you use an electrosurgical unit and an endoscope that is supplied with energy at the same time, the low-frequency currents from both power sources accumulate to form one single patient leakage current. If the electrosurgical unit belongs to type CF and the endoscope to type BF, the less favorable limit values for type BF apply to the overall patient leakage current. You should therefore only use power sources with the same type classification, ideally those of type CF.

5.8 Risk of burns to the patient

Flash-overs, activation in the trocar or without visual control of the distal end of the instrument can cause unwanted thermal injuries to the patient. Therefore, check that the distal end of the instrument protrudes at least 10 mm out of the trocar before each activation.

5.9 Danger of burns

Argon plasma can strongly heat electrically conductive materials, e.g. metal. To avoid burns, never direct argon plasma towards electrically conductive materials.

Risk of burns! Do not touch the instrument tip during activation (e.g. to clean the electrode).

The instrument tip may remain hot for a while after use. Do not touch the instrument tip and do not set it down on the patient or on any flammable materials.


5.10 Damage to the shaft tube

Do not bend the shaft tube in the case of applicators with a rotary wheel and moveable electrode.

In the case of the applicator with REF no. 20132-256 (with a flexible shaft tube): Only bend the shaft tube by hand, do not use any tools.

5.11 Modification of the product by the user

Erbe Elektromedizin expressly warns against modifying the product. Any modification exempts Erbe Elektromedizin from any and all liability.

6 How to use

6.1 Checking the expiration date

The expiration date is indicated on the packaging. Do not use this product if the expiration date has been exceeded!

6.2 Checking and opening the packaging

Only use this product if the packaging is unopened and undamaged. If the packaging is opened or damaged, the product is no longer sterile and must not be used.

Use aseptic techniques when removing the product from the packaging!

6.3 Checking the product




Check the product for damage before use, including all insulation.

If damaged, do not use this product!

6.4 Connecting the product and selecting unit settings

IMPORTANT! The APC unit is equipped with different APC socket. Depending on the type of socket, you may or may not need an adapter to connect the applicator.


To determine your type of socket please read the chapter "Required accessories".

Connection to an APC 2 with APC socket :

IMPORTANT! This type of socket is available in two variants, with and without finger switch activation. For the applicator, you will need the variant with finger switch activation, which can be recognized by the large metal ring at the left contact of the socket.

1. Switch on the VIO electrosurgical unit and select the basic program.
2. Insert the APC connector (1) into the APC socket (2).


3. For the connected applicator, a standard setting for mode, power, effect, and gas flow is automatically made on the VIO electrosurgical unit. You can change this standard setting within a certain range as required.

Connection to an APC 2 with APC socket :

IMPORTANT! The adapter with REF no. 20132-249 must be connected to the APC socket!

1. Switch on the VIO electrosurgical unit and select the basic program.
2. Insert the APC connector (1) into the adapter (2).

3. For the connected applicator, a standard setting for mode, power, effect, and gas flow is automatically made on the VIO electrosurgical unit. You can change this standard setting within a certain range as required.

Connection to an APC 3:

1. Switch on the VIO electrosurgical unit and pull up the main menu.
2. Insert the APC connector into the APC socket on the APC 3.
3. For the connected applicator, a standard setting for mode, effect level and gas flow is automatically made on the VIO electrosurgical unit. You can change this standard setting within a certain range as required.

6.5 In the case of endoscopic application: Selecting the trocar

Select a trocar for instruments with a diameter of 5 mm.

6.6 Using the product

During use, check the product for damage, including all insulation.

Press the yellow button to activate cutting mode. Press the blue button to activate coagulation mode. Make sure that the electro-surgical unit responds correctly.




Using the light gray ReMode button, you can, in the case of a program with two setting options *a* and *b*, toggle between both setting options.

In order to move the electrode longitudinally, turn the rotary wheel, or alternatively, move the shaft tube back and forth.

Spatula electrodes can be turned freely in any position. To do so, turn the shaft tube.



In the case of laparoscopic application: Retract the electrode completely when introducing the applicator or removing it from the trocar. You must feel the rotary wheel clicking into place. Only introduce the applicator into the trocar under visual control.

REF no. 20132-256: This applicator has a flexible shaft tube (**WARNING!** Only bend by hand, do not use any tools) and only one blue button. The blue button activates the coagulation mode. Make sure that the electrosurgical unit responds correctly.

7 Required accessories

IMPORTANT! The APC unit is equipped with different APC socket. Depending on the type of socket, you may or may not need an adapter to connect the applicator.

Socket type	Accessories for connecting the applicator
	No further accessories required.
	Adapter REF no. 20132-249

8 Disposal

Dispose of the product, packaging material and accessories (if available) in accordance with the national guidelines and legislation applicable in each case.

9 Symbols

Symbol	Explanation	Symbol	Explanation
	Consult instructions for use		Caution
	Catalogue number		Batch code
	Manufacturer		Use-by date ¹⁾
	Keep away from sunlight		Keep dry
	Quantity (x)		Sterilized using ethylene oxide
	Do not use if package is damaged and consult instructions for use		Do not re-use
	Rigid shaft tube with needle		Rigid shaft tube with scoop
	Flexible shaft tube		European conformity marking

1.) The date next to the icon is in the format Year-Month-Day.



NOTES ON USE

USA

Rigid Argon Plasma Coagulation Applicators

20132-250, 20132-251, 20132-252, 20132-253, 20132-254,
20132-255, 20132-256

STERILE EO



Table of Contents

IMPORTANT!	33
1 Product description.....	33
2 Intended Use	34
3 Normal Use	34
4 Maximum Electrical Capacity	34
5 Safety Instructions	35
6 How to use.....	39
6.1 Check the expiration date	39
6.2 Control and open packaging	39
6.3 Check the product	39
6.4 Connect the product and select Unit settings	40
6.5 In the case of endoscopic application: Selecting the trocar	41
6.6 Using the product.....	41
7 Disposal	42
8 Symbols	42


IMPORTANT!

Please read all information carefully.

These notes on use do not replace the user manual of the APC unit and Electrosurgical unit used! Read the user manuals of the APC unit and Electrosurgical unit and ask Erbe or your distributor in case of doubt!

Caution! Federal law (USA) restricts this device to sale by or on the order of a physician!

1 Product description

- 1 Yellow button: Activates CUT
- 2 Blue button: Activates COAG
- 3 Light gray ReMode key: Toggles between two setting options of a program (only possible with the following unit combinations: APC 2 with VIO 300 D, APC 3 with VIO 3)
- 4 Rotary wheel: for freely adjustable movement of the HF electrode
- 5 Shaft tube: for turning spatula electrodes
- 6 APC connector: for connection to the APC unit


2 Intended Use

20132-256: The Rigid Argon Plasma Coagulation Applicator is intended for monopolar coagulation and ablation of tissue using argon plasma. It is used in open surgical procedures.

20132-250, 20132-251, 20132-252, 20132-253, 20132-254 and 20132-255: The Rigid Argon Plasma Coagulation Applicators are intended for monopolar cutting and coagulation of tissue with and without argon gas or for coagulation and ablation using argon plasma. They are used in open and in laparoscopic surgery.

- used in open surgery 20132-250, 20132-251, 20132-252, 20132-253.
- used in laparoscopic surgery 20132-254, 20132-255.


3 Normal Use

The argon plasma coagulation applicator is intended for use with an APC 2 / APC 3 in combination with an Erbe VIO electrosurgical unit.

The APC 2 must be equipped with an APC socket with finger switch activation. The APC 3 is already equipped with a suitable socket ex works.


The product is intended for single use.

4 Maximum Electrical Capacity

The maximum electrical capacity of these products is:

- 4300 V_p

5 Safety Instructions

WARNING!

5.1 Sterile product for single use

Erbe recommends that the product should not be processed after use. Processing can affect the material properties and/or function of the product such that use in accordance with the intended use is no longer possible. Also, it cannot be ensured that the product will be sterile after processing. Risks include:

- Components/connections that are damaged by high temperatures (e.g. in the autoclave).
- Components (e.g., filter, extremely small lumen, etc.) that may prevent internal cleaning/disinfection.

5.2 User requirements

This product may only be used by medical personnel who have been trained per these instructions.

Do not use the product if you are unfamiliar with its operation in conjunction with the ESU and the Argon Plasma Coagulator.

5.3 General information about protection of the user, patient, and product

Do not use in the presence of combustible or explosive materials!

Never lay this product on the patient or in his/her direct vicinity.

In the case of a long activation time and/or activations which follow in quick succession (with no cooling time) the distal instrument tip can become extremely hot. This is not a material fault, rather a consequence of heavy-duty use.

Protect the product from all physical damage! Do not throw! Do not apply force of any kind!

Do not kink the tube or wrap it round the product to which it belongs!

If damaged, do not use this product!

Checks other than those described in these notes on use do not need to be carried out.



5.4 Danger of fire in bodily cavities

Before using the APC always ensure that there are no combustible or fire-supporting gases in the lumen of the organ involved.

5.5 Gas embolisms, gas emphysema


To prevent gas embolisms, the argon flow rate must not be set high enough that argon is blown into open vessels. To prevent gas embolisms or emphysema, never point the distal end of APC applicators/APC probes directly toward open vessels or press against tissue.

5.6 Danger of unintentional damage to nerve structures

Bear in mind that during an APC application in the direct vicinity of nerve structures stimulation may be induced by electric currents.

In the case of nerves which require protection it is advisable to keep the APC activation time short and adapt maximum power output to the required indication.

5.7 Risks in the case of endoscopic application

Risk of capacitive coupling: HF current can be transferred from the activated instrument to the trocar, even if the two are not touching (known as "capacitive coupling"). Risk of burns to the patient and/or operating theater staff. Do not use hybrid trocars (e.g. a combination of metal and plastic). Only use trocars designed for use with electrosurgical devices.

Unwanted stimulation of nerves and muscles due to low-frequency currents. Dangerous summation of patient leakage currents from several power sources: It is not possible to entirely avoid low-frequency currents in electrosurgery, but they can be limited by applying technical measures to equipment and accessories. Limit types show how wide-ranging these measures are, e.g., type B, BF or CF as per EN 60601. Type CF offers the highest level of protection against low-frequency currents. However, please note the following: If you use an electrosurgical unit and an endoscope that is supplied with energy at the same time, the low-frequency currents from both power sources accumulate to form one single patient leakage current. If the electrosurgical unit belongs to type CF and the endoscope to type BF, the less favorable limit values for type BF apply to the overall patient leakage current. You should therefore only use power sources with the same type classification, ideally those of type CF.

5.8 Risk of burns to the patient

Flash-overs, activation in the trocar or without visual control of the distal end of the instrument can cause unwanted thermal injuries to the patient. Therefore, check that the distal end of the instrument protrudes at least 10 mm out of the trocar before each activation.

5.9 Danger of burns

Argon plasma can strongly heat electrically conductive materials, e.g. metal. To avoid burns, never direct argon plasma towards electrically conductive materials.

Risk of burns! Do not touch the instrument tip during activation (e.g. to clean the electrode).

The instrument tip may remain hot for a while after use. Do not touch the instrument tip and do not set it down on the patient or on any flammable materials.

5.10 Damage to the shaft tube

Do not bend the shaft tube in the case of applicators with a rotary wheel and moveable electrode.

In the case of the applicator with REF no. 20132-256 (with a flexible shaft tube): Only bend the shaft tube by hand, do not use any tools.

5.11 Modification of the product by the user

Erbe Elektromedizin expressly warns against modifying the product. Any modification exempts Erbe Elektromedizin from any and all liability.

6 How to use


6.1 Check the expiration date

An expiration date is printed on the packaging of the product. Do not use any product beyond its expiration date!

6.2 Control and open packaging

Use this product only if the packaging has not been opened or damaged. If the packaging is open or damaged, the product is no longer sterile and must not be used. In this case, throw the product away!




Use aseptic techniques when removing the product from the packaging!

6.3 Check the product

Check the product for damage before use, including all insulation.

If damaged, do not use this product!


6.4 Connect the product and select Unit settings

Connection to an APC 2 with APC socket :

IMPORTANT! This type of socket is available in two variants, with and without finger switch activation. For the applicator, you will need the variant with finger switch activation, which can be recognized by the large metal ring at the left contact of the socket.

1. Turn on the ESU (i.e., the VIO Model) and select the basic program.
2. Plug the APC connector (1) of the applicator into the APC socket (2).

3. For the connected applicator, a standard setting for mode, power, effect, and gas flow is automatically made on the VIO electrosurgical unit. You can change this standard setting within a certain range as required.

Connection to an APC 3:

1. Switch on the VIO electrosurgical unit and pull up the main menu.
2. Insert the APC connector into the APC socket on the APC 3.
3. For the connected applicator, a standard setting for mode, effect level and gas flow is automatically made on the VIO electrosurgical unit. You can change this standard setting within a certain range as required.

6.5 In the case of endoscopic application: Selecting the trocar

Select a trocar for instruments with a diameter of 5 mm.

6.6 Using the product

During use, check the product for damage, including all insulation.

Press the yellow button to activate cutting mode. Press the blue button to activate coagulation mode. Make sure that the electro-surgical unit responds correctly.




Using the light gray ReMode button, you can, in the case of a program with two setting options *a* and *b*, toggle between both setting options.

In order to move the electrode longitudinally, turn the rotary wheel, or alternatively, move the shaft tube back and forth.

Spatula electrodes can be turned freely in any position. To do so, turn the shaft tube.

In the case of laparoscopic application: Retract the electrode completely when introducing the applicator or removing it from the trocar. You must feel the rotary wheel clicking into place. Only introduce the applicator into the trocar under visual control.

REF no. 20132-256: This applicator has a flexible shaft tube (**WARNING!** Only bend by hand, do not use any tools) and only one blue button. The blue button activates the coagulation mode. Make sure that the electrosurgical unit responds correctly.

7 Disposal

Dispose of the product, packaging material and accessories (if available) in accordance with the national guidelines and regulation applicable in each case.

8 Symbols

ISO 15223-1 Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied

Ref. no.	Symbol	Title	Explanation
5.4.3		Consult instructions for use	Indicates the need for the user to consult the instructions for use.
5.1.6		Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
5.1.1		Manufacturer	Indicates the medical device manufacturer.
5.3.2		Keep away from sunlight	Indicates a medical device that needs protection from light sources.
5.3.4		Keep dry	Indicates a medical device that needs to be protected from moisture.


Ref. no.	Symbol	Title	Explanation
5.4.4		Caution	Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.
5.1.5		Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
5.1.4		Use-by date ¹⁾	Indicates the date after which the medical device is not to be used.
5.2.3		Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide.
5.2.8		Do not use if package is damaged and consult instructions for use	Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information.
5.4.2		Do not re-use	Indicates a medical device that is intended for one single use only.

1.) The date next to the icon is in the format Year-Month-Day.

MDD 93/42/EEC, Annex XII

Symbol	Title	Explanation
	European conformity marking	Confirmation from the manufacturer that the product meets the requirements of the applicable European guidelines.

Symbols not relating to standards

Symbol	Title	Explanation
	Quantity (x)	Quantity of the product. x pieces of the product are supplied in one packaging unit.
	Rigid shaft tube with needle	This instrument has a rigid, non-bendable instrument tip with a needle electrode.
	Rigid shaft tube with scoop	This instrument has a rigid, non-bendable instrument tip with a spatula electrode.
	Flexible shaft tube	This instrument has a bendable instrument tip.

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