
erbe
power your performance.

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- 
- 1 Kaehler GF, Sold MG, Post S, Fischer K, Enderle MD. Selective Tissue Elevation by Pressure Injection (STEP) Facilitates Endoscopic Mucosal Resection (EMR). *Surg Technol Int* 2007; 16: 107–12
- 10 Zhou PH, Schumacher B, Yao LQ, Xu MD, Nordmann T, Cai MY, Charton JP, Vieth M, Neuhaus H. Conventional vs. water-jet-assisted endoscopic submucosal dissection in early gastric cancer: a randomized controlled trial. *Endoscopy* 2014; 46: 836–842
- 12 Schumacher B, Charton JP, Nordmann T, Vieth M, Enderle MD, Neuhaus H. Endoscopic submucosal dissection of early gastric neoplasia with a water jet-assisted knife: a Western, single-center experience. *Gastrointest Endosc.* 2012 Jun; 75(6): 1166–74
- 18 Cai MY, Zhou PH, Yao LQ, Xu MD, Zhong YS, Li QL, Chen WF, Hu JW, Cui Z, Zhu BQ. Peroral endoscopic myotomy for idiopathic achalasia: randomized comparison of water-jet assisted versus conventional dissection technique. *Surg Endosc.* 2014 Apr; 28(4): 1158–65

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- 1 Kaehler GF, Sold MG, Post S, Fischer K, Enderle MD. Selective Tissue Elevation by Pressure Injection (STEP) Facilitates Endoscopic Mucosal Resection (EMR). *Surg Technol Int* 2007; 16: 107–12
- 2 Belle S, von Boscamp M, Sold M, Mack S, Pilz L, Ebert M, Kaehler G. Submucosal injection with waterjet improves endoscopic mucosal resection of colorectal adenoma – a randomised controlled clinical trial. *Scand J Gastroenterol* 2017; VOL. 52, NO. 2, 222–227
- 2a Belle S, Collet PH, Szyrach M, Ströbel P, Post S, Enderle MD, Kähler G. Selective tissue elevation by pressure for endoscopic mucosal resection of colorectal adenoma: first clinical trial. *Surgical Endoscopy* 2012; Volume 26

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

- 10 Zhou PH, Schumacher B, Yao LQ, Xu M-D, Nordmann T, Cai MY, Charton JP, Vieth M, Neuhaus H. Conventional vs. water-jet-assisted endoscopic submucosal dissection in early gastric cancer: a randomized controlled trial. *Endoscopy* 2014; 46: 836–842
- 11 Repici A, Hassan C, Pagano N, Rando G, Romeo F, Spaggiari P, Ron-calli M, Ferrara E, Malesci A. High efficacy of endoscopic submucosal dissection forectal laterally spreading tumors larger than 3 cm. *Gastrointest Endosc.* 2013 Jan; 77(1): 96–101
- 12 Schumacher B, Charton JP, Nordmann T, Vieth M, Enderle MD, Neuhaus H. Endoscopic submucosal dissection of early gastric neoplasia with a water jet-assisted knife: a Western, single-center experience. *Gastrointest Endosc* 2012, Jun; 75(6): 1166–74
- 13 Repici, A et al. A prospective, single center study of endoscopic sub-mucosal dissection of rectal LST lesion larger than 3cm by using an innovative concept of injecting and Cutting: The water-Jet Hybrid-Knife (ESD-H). *Gastrointest Endosc Vol 73, Issue 4, Supplement, Page AB156, April 2011*
- 14 Neuhaus, H. Endoscopic mucosal resection and endoscopic submucosal dissection in the West – too many concerns and caveats? *Endoscopy* 2010; 42: 859–861

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HybridKnife® is not just elevating, it also allows me to manage all the other steps of an ESD with the same instrument. In combination with the modes of the VIO® 3 generator, HybridKnife® is my preferred choice for the en-bloc resection of lesions in the GI tract.

Prof. Dr. med. Alessandro Repici
Klinischer Direktor
Prof. of Gastroenterology Humanitas
University, Mailand

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I started using the O-type Hybridknife® next to the I-type in the more difficult to reach areas and learned about the great versatility this knife has to offer. It really makes ESD a single-instrument procedure. It makes a full circumferential incision in a 'lift-and-cut' fashion, without having to remove the knife a single time, rather easy. This is true even in retroflex positions.


The knife can cut in forward positions as well as lateral cutting and dissecting while hooking tissue in awkward positions without having to fear deeper thermal injury. Bleedings can be managed with the knife with submucosal lifting at my disposal at any time. I hardly use a coagulation forceps anymore.


Dr. med. Arjun Koch
Department of Gastroenterology and Hepatology
Erasmus MC, Rotterdam

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

17 Li QL, Zhou PH. Perspective on peroral endoscopic myotomy for achalasia: Zhongshan experience. *Gut Liver.* 2015 Mar; 9(2): 152–8

18 Cai MY, Zhou PH, Yao LQ, Xu MD, Zhong YS, Li QL, Chen WF, Hu JW, Cui Z, Zhu BQ. Peroral endoscopic myotomy for idiopathic achalasia: randomized comparison of water-jet assisted versus conventional dissection technique. *Surg Endosc.* 2014 Apr; 28(4): 1158–65

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HybridKnife® is my favorite knife in endoscopic resection, such as ESD and POEM procedures for its efficient submucosal injection for creating the submucosal tunnel and accurate cutting and coagulation.

Prof. Dr. med. Ping-Hong Zhou
Head of Endoscopy Center
Fudan University Zhongshan Hospital, Shanghai

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Thanks to HybridKnife®, I can significantly reduce the time taken for doing an endoscopic surgical procedure especially when there is fibrosis and hypervascularity. The credit goes to the all in one capability of this wonderful instrument.

Dr. med. Amit Prabhaka Maydeo
Baldota Institute of Digestive Sciences
Global Hospital, Mumbai

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Erbe Elektromedizin GmbH