
erbe
power your performance.

ENDO CUT® I

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Arvanitakis M, Devière J. Endoscopic Retrograde Cholangiopancreatography. *Endoscopy* 2004; 36:855-859.

ASGE guideline: the role of ERCP in disease of the biliary tract and the pancreas. *Gastrointest Endosc* 2005, 62: 1-8.

Baillie J. Complications of ERCP. ERCP and Its Applications, edited by I.M. Jacobson, Lippincott Raven Publishers, Philadelphia, 1998.

Binmoeller K, Seifert H, Gerke H, Seitz U, Portis M, Soehendra N. Papillary roof incision using the Erlangen-type precut papillotome to achieve selective bile duct cannulation. *Gastrointest Endosc* 1996; 44:689-695.

Bruins Slot W, Schoeman MN, Disario JA, Wolters F, Tytgat GN, Huibregtse K. Needle-knife sphincterotomy as a precut procedure: a retrospective evaluation of efficacy and complications. *Endoscopy* 1996; 28:334-339.

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Classen M, Born P. Endoscopic sphincterotomy. In: Tytgat GNJ, Classen M, Wayne JD, Nakazawa S, editors. Practice of therapeutic endoscopy. 3rd ed. Philadelphia: Saunders, 2000:129-46.

Costamagna G. State of the Art Review – Therapeutic biliary endoscopy. *Endoscopy* 2000; 32:209-216.

Cotton PB, Lehman G, Vennes J, Geenen JE, Russell RC, Meyers WC, et al. Endoscopic sphincterotomy complications and their management: an attempt at consensus. *Gastrointest. Endosc.* 1991; 37:383-93.

Cotton PB. Precut sphincterotomy: a risky technique for the experts only. *Gastrointest Endosc* 1989; 35:578-579.

DiSario JA, Freeman ML, Bjorkman DJ, Macmathuna P, Petersen BT, Jaffe PE, Morales TG, Hixson LJ, Sherman S, Lehman GA, Jamal MM, Al-Kawas FH, Khandelwal M, Moore JP, Derfus GA, Jamidar PA, Ramirez FC, Ryan ME, Woods KL, Carr-Locke DL, Alder SC. Endoscopic balloon dilation compared with sphincterotomy for extraction of bile duct stones. *Gastroenterology* 2004; 127:1291-1291.

Elfant AB, Bourke MJ, Alhalel R, Kortan PP, Haber GB. A prospective study of the safety of endoscopic therapy for choledocholithiasis in an outpatient population. *Am J Gastroenterol* 1996; 91:1499-1502.

Ell C, Rabenstein T, Ruppert T, Forster P, Hahn EG, Demling L. 20 years of endoscopic papillotomy: analysis of 2752 patients in Erlangen Hospital. *Dtsch Med Wochenschr* 1995; 120:163-7.

Elta GH, Barnett JL, Wille RT, Brown KA, Chey WD, Scheiman JM. Pure cut electrocautery current for sphincterotomy causes less post-procedure pancreatitis than blended current. *Gastrointest Endosc.* 1998 Feb;47(2):149-53.

Fogel E, McHenry L, Sherman S, Watkins J, Lehman G. Therapeutic Biliary Endoscopy. *Endoscopy* 2005; 37:139-145.

Foutch PG. A prospective assessment of results for needle-knife papillotomy and standard endoscopic sphincterotomy. *Gastrointest Endosc* 1995; 41:25-32.

Fuji T, Amano H, Ohmura R, Akiyama T, Aibe T, Takemoto T. Endoscopic pancreatic sphincterotomy: technique and evaluation. *Endoscopy* 1989; 21:27-30.

Fujita N, Maguchi H, Komatsu Y, Yasuda I, Hasebe O, Igarashi Y, Murakami A, Mukai H, Fujii T, Yamao K, Maeshiro K; JESED Study Group. Endoscopic sphincterotomy and endoscopic papillary balloon dilatation for bile duct stones: a prospective randomized controlled multicenter trial. *Gastrointest Endosc* 2003; 57:151-155.

- Freemann ML, Guda NM. ERCP cannulation: a review of reported techniques. *Gastrointest Endosc* 2005; 61:112-125.
- Freemann ML. Adverse outcomes of endoscopic retrograde cholangiopancreatography: avoidance and management. *Gastrointest Endosc Clin N Am* 2003; 13:775-98, xi.
- Freeman ML, Guda NM. Prevention of Post-ERCP pancreatitis: a comprehensive review. *Gastrointest Endosc* 2004; 59:845-864.
- Freeman ML. VI. Adverse events and success of ERCP – Adverse outcomes of ERCP. *Gastrointest Endosc* 2002; 56:273-282.
- Freeman ML, DiSario JA, Nelson DB, et al. Risk factors for post-ERCP pancreatitis: a prospective, multicenter study. *Gastrointest Endosc* 2001; 54:425-34.
- Freeman M. Complications of endoscopic biliary sphincterotomy: a review. *Endoscopy* 1997; 29:288-97.
- Freeman ML, Nelson DB, Sherman S, Haber GB, Herman ME, Dorsher PJ, et al. Complications of endoscopic biliary sphincterotomy. *N Engl J Med* 1996; 335:909-918.
- Goff JS. Long-term experience with the transpancreatic sphincter pre-cut approach to biliary sphincterotomy. *Gastrointest Endosc* 1999; 50:642-645.
- Gorelick A, Cannon M, Barnett J, Chey W, Scheiman J, Elta G. First cut, then blend: An electrocautery technique affecting bleeding at sphincterotomy. *Endoscopy* 2001; 33:976-980.
- Guitron A, Adalid R, Barinagarrementeria R, Gutierrez JA, Mena G, Nares J. Pre-cut sphincterotomy: efficacy and complications. *Rev Gastroenterol Mex* 1998; 63:148-52.
- Guitron A, Adalid R. Safety and efficacy of sphincterotomy using guided and pre-cut techniques. A 5-year experience. *Rev Gastroenterol Mex* 1996; 61:342-1347.
- Häfner M and Schöfl, R. Diagnostic retrograd cholangiopancreatograph
- Hashiba K, D'Assunção MA, Armellini S, Hassegawa RT, Cappellanes CA, Moribe D. Endoscopic Suprapapillary Blunt Dissection of the Distal Common Bile Duct in Cases of Difficult Cannulation: A Pilot Series. *Endoscopy* 2004;36: 717-321.
- Hochberger J. Difficult Cannulation. In: Ginsberg G, Kochman M, Norton I, Gostout CJ, eds. *Clinical Gastrointestinal Endoscopy. Volume 1: Elsevier Saunders, USA, 2005:617-24.*
- Jacobson I. ERCP and its Applications. 1998; Lippincott – Raven Publishers, Philadelphia, New York.
- Kaffes AJ, Sriram PV, Rao GV, Santosh D, Reddy DN. Early institution of pre-cutting for difficult biliary cannulation: a prospective study comparing conventional vs. a modified technique. *Gastrointest Endosc*. 2005; 62:669-674.
- Kasmin FE, Cohen D, Batra S, Cohen SA, Siegel JH. Needle-knife sphincterotomy in a tertiary referral center: efficacy and complications. *Gastrointest Endosc* 1996; 44:48-53.
- Kohler A, Maier M, Benz C, Martin WR, Farin G, Riemann JF. A new HF current generator with automatically controlled system (Endocut mode) for endoscopic sphincterotomy: preliminary experience. *Endoscopy* 1998; 30:351-355.
- Lella F, Bagnolo, Colombo E, Umberto Bonassi, MD, et al. A simple way of avoiding post-ERCP pancreatitis. *Gastrointest Endosc*. 2004; 59:830-834.
- Levy MJ, Geenen JE. Endoscopic management of benign pancreatic diseases. In: Tytgat GN, Classen M, Waze JD, Nakazawa J, editors. *Practive of tehrapetic endoscopy*. 2nd ed. London: Saunders, 2000:177-197.
- Macintosh DG, Love J, Abraham NS. Endoscopic sphincterotomy by using pure-cut electrosurgical current and the risk of post-ERCP pancreatitis: a prospective randomized trial. *Gastrointest. Endosc*. 2004; 60:551-556.
- Martin DF, England R, Martin O. The safety sphincterotome: the device, the technique, and preliminary results. *Endoscopy* 1998; 30:375-378.
- Masci E, Toti G, Mariani A, et al. Complications of diagnostic and therapeutic ERCP: A prospective multicenter study. *Am J Gastroenterol* 2001; 96:417-423.
- Maguchi H, Takahashi K, Kanatuma A. Pancreatic and bile duct stricture – Endoscopic management of biliary and pancreatic ducts strictures. *Digest Endosc* 2004; 16 (Suppl.):34-38.
- Maydeo A, D. Borkar D. Techniques of selective cannulation and sphincterotomy. *Endoscopy* 2003; 35:19-23.
- Mavrogiannis C, Liatsos C, Romanos A, Petoumenos C, Nakos A, Karvountzis G. Needle-knife fistulotomy versus needle-knife precut papillotomy for the treatment of common bile duct stones. *Gastrointest Endosc* 1999; 50:334-9.
- Mergener K, Kozarek RA. Therapeutic pancreatic endoscopy. *Endoscopy* 2005; 37:201-207.
- Neuhaus H. Therapeutic pancreatic endoscopy. *Endoscopy* 2004; 36: 8-16.
- Norton ID, Petersen BT, Bosco J, Nelson DB, Meier PB, Baron TH, Lange SM, Gostout CJ, Loeb DS, Levy MJ, Wiersema MJ, Pochron N. A randomized trial of endoscopic biliary sphincterotomy using pure-cut versus combined cut and coagulation waveforms. *Clin Gastroenterol Hepatol*. 2005; 3:1029-1033.
- Park SH, Kim HJ, Park DH, Kim JH, Lee JK, Lee SH, Chung IK, Kim HS, Kim SJ. Pre-cut papillotomy with a new papillotome. *Gastrointest Endosc* 2005; 62: 588-591.
- Kalloor AN, Thuluvath PJ, Pasricha PJ. Treatment of high-risk patients with symptomatic cholelithiasis by endoscopic gallbladder stenting *Gastrointest Endosc*. 1994 Sep-Oct;40(5):608-10.
- Perini RF, Sadurski R, Cotton PB, Patel RS, Hawes RH, Cunningham JT. Post-sphincterotomy bleeding after the introduction of microprocessor-controlled electrosurgery: does the new technology make the difference? *Gastrointest Endosc* 2005; 61:53-57.
- Rabenstein T, Schneider HT, Nicklas M, Ruppert T, Katalinic A, Hahn EG, et al. Impact of skill and experience of the endoscopist on the outcome of endoscopic sphincterotomy techniques. *Gastrointest Endosc* 1999; 50:628-36.
- Rabenstein T, Ruppert T, Schneider HT, Hahn EG, Ell C. Benefits and risks of needle-knife papillotomy. *Gastrointest Endosc* 1997; 46:207-11.
- Riemann JF. [Endoscopic treatment of gallstones: new technics] *Z Gesamte Inn Med*. 1989 Oct 15;44(20):605-9.
- Rollhauser C, Johnson M, Al-Kawas FH. Needle-knife papillotomy: a helpful and safe adjunct to endoscopic retrograde cholangiopancreatography in a selected population. *Endoscopy* 1998; 30:691-696.
- Sherman S, Ruffolo T, Hawes RH, Lehman GA. Complications of endoscopic sphincterotomy. *Gastroenterology* 1991; 101: 1068-1075.
- Silverstein F.E., Tytgat G.N.J. *Praxis der Gastroenterologischen Endoskopie; Atlas und Lehrbuch*, 1999; Georg Thieme Verlag, Stuttgart
- Soehendra N, Binmoeller KF, Seifert H, Schreiber HM. *Therapeutic Endoscopy. Color Atlas of Operative Techniques for the Gastrointestinal Tract*. 2005, Georg Thieme Verlag, Stuttgart.
- Sriram PV, Rao GV, Nageshwar Reddy D. The precut: when, where and how? A review. *Endoscopy* 2003; 35: 24-30.
- Stefanidis G, Karamanolis G, Viazis N, Sgouros S, Papadopoulos E, Ntatsakis K, Mantides A, Nastos H. A comparative study of postendoscopic sphincterotomy complications with various types of electrosurgical current in patients with choledocholithiasis. *Gastrointest Endosc* 2003; 57: 192-197.
- Sugiyama M, Suzuki Y, Abe N, Masaki T, Mori T, Atomi Y. Pancreas and Biliary – Endoscopic retreatment of recurrent choledocholithiasis after sphincterotomy. *Gut* 2004; 53:1856-1859.
- Tang SJ, Haber GB, Kortan P, Zanati S, Cirocco M, Ennis M, Elfant A, Scheider T, Ter H, J. Dorais J. Precut papillotomy versus persistence in difficult biliary cannulation: A prospective randomized trial. *Endoscopy* 2005; 37:58-65.

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Erbe Elektromedizin GmbH