

FOR RESECTION AND VAPORIZATION IN SALINE SOLUTION:
BIPOLAR CUT++.

1 BIPOLAR RESECTION

Guide / progs. MF receptacle

Mode BIPOLAR CUT++	Mode BIPOLAR SOFT++
Effect 5	Effect 6
Upmax:390Vp	Upmax:155Vp




BIPOLAR CUT++: THE FIRST PLUS FOR CUTTING AND THE SECOND PLUS FOR HEMOSTASIS.

BIPOLAR CUT++ offers a new cutting quality for resection of the prostate, bladder and uterus.

And these are the advantages:

- ✦ Rapid generation of arcs, immediate cutting (without needing to exert mechanical pressure on the tissue)
- ✦ Arcing intensity is regulated to ensure homogeneous development of arcs (at low energy level inputs)
 - consistent vaporization effects
 - safe for patients, homogeneous hemostasis
 - reduces operative costs (as fewer snares are required)
- ✦ User-friendly: plug and play
- ✦ BIPOLAR CUT++ is compatible with all major bipolar resectoscopes

With AUTO CUT, HIGH CUT and DRY CUT the VIO system also provides a range of tried and tested monopolar resection modes.



Upgrade your VIO 300 D with the bipolar resection adapter and appropriate software



BIPOLAR CUT++ with automatic output regulation for bipolar resection in saline solution.



Rapid generation of arcs.



Homogeneous effects for the vaporization of tissue structures.

ERBE Elektromedizin GmbH
Waldhoernlestrasse 17
72072 Tuebingen
Germany
Phone +49 7071 755-0
Fax +49 7071 755-179
sales@erbe-med.de
www.erbe-med.com



ERBE

Perfection for Life